



35 Industrial Parkway, Somerville, NJ. 08876 U.S.A.

Tel: 908-722-5000 Fax: 908-722-5024

www.biosearch.com

HOSPITAL PRICE LIST (2006)

ANORECTAL BIOFEEDBACK MONITORS FOR FECAL INCONTINENCE AND FECAL CONSTIPATION


The Biosearch® Anorectal Biofeedback Devices are hand held, home training aids for patients who suffer from either fecal incontinence (Model 5) or fecal constipation (Model 10). Both devices are used for anal sphincter muscle re-education and consist of two main components; the anorectal probe and the monitor. The silicone anorectal balloon probe consists of two balloons - the distal balloon which is filled with air to simulate the urge to defecate and the proximal balloon which is a pneumatic sensor for anal manometry. The rechargeable hand-held monitor provides a visual representation of the patients anal resting and squeezing activity by converting the pneumatic signal from the probe into a DC signal which illuminates the LED bar graph on the monitor. What makes our device so unusual is it's very low cost and it's simplicity, thus allowing most patients to afford and operate our device in the privacy of their own homes.

{1 unit per box}

<u>Reorder No.</u>	<u>Description</u>	<u>Price</u>
18-1005	(Biofeedback monitor #5 & 1 probe) (Fecal Incontinence)	\$650.00
18-1010	(Biofeedback monitor #10 & 1 probe) (Fecal Constipation)	\$650.00
18-1000	(Replacement probe)	\$125.00

Hand Held
Anorectal Biofeedback

*Strip Chart
Recorder*




*System 5
For Fecal
Incontinence*



*System 10
For Chronic
Constipation*



*Anorectal
Probe*



- High tech ergonomically designed monitor is simple to use and provides accurate results.
- Illuminated color LEDs indicate patient biofeedback response.
- Medical grade silicone anorectal probe is anatomically designed for patient comfort and biocompatibility.
- Optional Strip Chart Recorder supplies a permanent record of training sessions.

Anorectal Biofeedback

HAND HELD BIOFEEDBACK MONITORS

- Sensitivity adjustment adaptable to patients "squeeze" ability.
- Rechargeable battery and recharger supply long lasting power.
- Color LED visual and/or audible signals assist sphincter training.
- Affordably priced for both the physician and patient.
- Compact design with discrete carrying case for easy transport.
- Approximate dimensions: 5 $\frac{5}{8}$ "H x 3 $\frac{1}{4}$ "W x $\frac{7}{8}$ "D.
- Approximate weight = 8.8 ounces.

STRIP CHART RECORDER (OPTIONAL)

- Inkless, dry writing on pressure sensitive chart paper.
- Access window slides open so notations can be easily made on chart paper.
- 6 cm (2.3) chart width.
- Single channel.

ANORECTAL PROBE

- Anatomically designed probe responsive to sphincter contraction and relaxation.
- Single patient use avoids contamination.
- Easy to clean with soap and water.
- Replacement probes are available.



FOR MORE INFORMATION ON THIS PRODUCT,
CONTACT YOUR LOCAL BIOSEARCH DEALER.

Description	Product Number
Anorectal Biofeedback System 5 Kit (Fecal Incontinence) includes: 1 - Biofeedback monitor, model 5 1 - Anorectal probe	18-1005
Anorectal Biofeedback System 10 Kit (Constipation) includes: 1 - Biofeedback monitor, model 10 1 - Anorectal probe	18-1010
Anorectal Probe (replacement)	18-1000
Strip Chart Recorder (optional)	18-1900
Chart Paper (replacement)	18-1901

biosearch[®]

The Medical Innovator[™]

Biosearch Medical Products, Inc.
P.O. Box 1700, 35 Industrial Parkway
Somerville, NJ 08876-1276 U.S.A.
(800) 326-5976 (908) 722-5000 Fax (908) 722-0132

Visit our Web site at <http://www.mc2usa.com/biosrch/>