

biosearch[®]

MEDICAL PRODUCTS INC.

35 Industrial Parkway, Somerville, NJ. 08876 U.S.A.

Tel: 908-722-5000 Fax: 908-722-5024

www.biosearch.com

(2007)

NASOGASTRIC FEED CATHETERS

Biosearch Medical Products, Inc. has been a leading manufacturer of polyurethane Nasogastric Feeding catheters for 30 years. Our catheters were the first to be made of Polyurethane which offers larger lumens (Fluid pathway) and smaller exterior diameters together with softness and high strength. A very important addition feature is that these Nasogastric Feeding catheters are coated with the Hydromer[®] hydrophilic lubricant which becomes slippery when made wet. This Hydromer coating is applied to the weight assemble to enhance ease of insertion and to the lumens to aid in Stylet removal.

We prefer to sell these to our clients in a “private label” configuration, with your name on the labels. Offer in either sterile (EO) or non-sterile packaging, sterile versions will require a minimum order of 1000 units (per delivery). Sold in bulk or in boxes of 10 units each.



{Price per unit}

<u>Reorder No.</u>	<u>Description</u>	<u>Price non-sterile</u>	<u>Price sterile</u>
BMP-7165	8ch/fr (2,7mm), 43 inches/109 cm long	\$Please	Call
BMP-7143	10ch/fr (3,3mm), 43 inches/109 cm long	\$Please	Call
BMP-7159	12ch/fr (4,0mm), 43 inches/109 cm long	\$Please	Call